

How to run software?

e10001/GG1

1. Run Readout

Type ./godaq

code location: /opt/lucid/daq/10.1-008/

run readout: ./godaq

2. Run Scalers display

code location: ~/VMUSB/Scalers

run: goscaler

Note: You can use "goscaler" from anywhere

3. Run Switcher

"goswitch" from anywhere

4. Run elog

code location: /user/e10001/elogsync

Type commands:

elogServer

-run this only ONCE. Running it on two different machines will kill the server

Goelog

Start elogclient; can be run multiple times.

5. Power supply printout for HiRA Si and CsI

telnet 35.9.56.159 1527

user name : admin

password: admin

Use tab to navigate the menus. Go to Main -> Channels, then Groups->Group 05.

DON'T CHANGE ANYTHING IF YOU DON'T KNOW WHAT YOU'RE DOING!!!

Printout the voltage log on linux

Applications→Accessories→Take screens shot

“Grab current window”, wait for a second

Click “Take screenshot” and immediately click the HV display window.

6. Online Data analysis with SpecTcl

code location: e10001/VMUSB/VMUSBSpecTcl

run spectcl : ./SpecTcl

files in /e10001/VMUSB/VMUSBSpectcl/

description	def-file [def-files]	win-file [win-files]
HiRA Si		
Summary of raw energy distribution for each Si	def-files/asic_energy_summary.tcl	win-files/asic_energy_summary.win
Calibrated energy summary		
Summary of time distribution for each Si	def-files/asic_time_summary.tcl	win-files/asic_time_summary.win
HiRA CsI		
Summary of raw energy distribution for each CsI	def-files/csi_summary.tcl	win-files/csi.win
CAESAR		
Summary	def-files/caesarsummary.tcl	win-files/CEsummary.win

For offline analysis:

Data Source → File → **click ringbuffer (very important or the program crashes)**

Runs are in stagearea/complete

7. HiRA Si Control Software

code location: ~/VMUSB/ASIC_control_E

run : ./CHIP

then

ONLY if the VME crate has been reset, either by software or by powering it off, load the XLM configurations a.k.a. the bit file. **DON'T DO THIS IF THE VME CRATE HAS NOT BEEN DISTURBED:**

- 1) On the File menu, select XLM configure
- 2) In the frame that opens at the right of the program, you will have to load a configuration into each XLM.
- 3) Set the "crate" slider to 0, and the "type" slider to XLMXXV
- 4) Set the "slot" slider to 3
- 5) push "Pick Load File", then select "Browse". You cannot just type the filename.
- 6) Select "xlmxxv_rev518.bit"
- 7) The bit file will load into the XLM, which will take some time (10 s to 1 m)
- 8) If the program crashes, chances are the XLM has locked up and you have to turn the VME crate off and back on again, and start over.

- 9) When it has loaded, check the messages in the terminal window. They will tell you whether the XLM correctly communicated with the motherboards. Make a note if you saw any errors here but you can continue regardless if there are only a couple of lines of errors.
- 10) Repeat steps 3 – 9 for the other two XLMs (slots 4 and 5).
- 11) From the File menu, select Load.
- 12) Select Browse. Again, you cannot just type in the filename
- 13) Browse to setupfiles, and select: MB123_0205.setup

DON'T CHANGE ANYTHING IF YOU DON'T KNOW WHAT YOU'RE DOING!!!

8. Load CsI

code location: ~/VMUSB/CSI_disc

run: wish pico.tcl

This will not work if CAMAC crate is off.

(Re)load gains and thresholds by pressing F8.