

Pulser Calibration

1. Load 1d_ene_all_RAW.tcl for data sorting
2. Save the spectra in asc format
3. Go to /projects/hira/lobastov/Pulser/
4. Type "PulserCal" to start program
5. Go to "File", click "open spectra asc file" then load the asc file.
6. Click "all spectra" to start automatic fitting
7. After fitting, save the vdef.files. Go to "File" and click "save offset/slope" and "save peak position".
8. To close the program, go to "File" and click "Exit"

Remark: Change PulserCal.cc for tower configuration and PulserCal.cfg for parameters and ranges if necessary.